NOTES

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USE CASES

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Post Processing WS

This system is intended to be used by health care professionals for the purpose of obtaining PVE corrected PET images and/or object values by using advanced post processing procedures.

Reporter

This actor is intended to be used for the purpose of :

- tracing the Post Processing Procedure procedures applied (including inputs used and outputs produced).

- caching the images involved in the process.

This actor will be implemented in the Rasna REPORTER environment.

The REPORTER is a multi-platform environment that provides full support for DICOM Structured Reporting.

It provides moreover tools to create an unambiguous semantic documentation for diagnosis, including report context, observer, previous reports, links with the medical study, images and other information.

The features that should be of interest for the project are :

- The interfaces for the import and the export of Structured Reports
- The image caching
- The Structured Reports storage

The available external interfaces (services) for the data exchange are:

- DICOM Storage Service Class. This allows to send and receive SRs and related images from and to external DICOM devices.
- HL7 export (ORU Messages) [Not applicable the project purpose]
- Web based services. Through these services external systems may import and export DICOM Structured report by using XML objects

SCENARIO PRINCIPALE

Descrizione


The PP WS processes PET/SPECT and MRI images in order to obtain the PVE corrected PET images and/or object values (Process Images) .

The Post Processing Procedure is traced and the results recorded (Trace PP Procedure)

Diagramma

Created/modified by: Giorgio Cangioli
Creation Date: 04/09/2002

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This document contains proprietary and confidential information of Rasna Imaging Systems. No part of this publication may be reproduced, transmitted, transcribed, recorded, or translated into any language in any form or by any means, electronic, magnetic, optical, chemical, physical or otherwise. Rasna Imaging Systems reserves the right to correct this publication and to make changes to its contents without the obligation of informing anyone of these revisions or changes. ©2002, Rasna Imaging Systems (09/02)

SCENARIO "PROCESS IMAGES"


Descrizione

The PP WS processes PET/SPECT and MRI images in order to obtain the PVE corrected PET images and/or object values .

This is obtained applying Co-Registration and Segmentation procedures before the PVE correction algorithm.

The PP WS makes use of the packages: Co-Registration, Segmentation and PVE Correction

Diagramma

SCENARIO "TRACE POST PROCESSING PROCEDURES"


Descrizione

The Post Processing Procedures applied are traced and the results recorded by using DICOM Structured Reports.

The "tracing" reports and the processing results are stored into the Reporter Repository.

Both the PP WS and the Reporter make use of the packages: Co-Registration, Segmentation and PVE Correction

Diagramma


CO-REGISTRATION

Descrizione

The PP WS gets the suitable inputs, performs the Co-Registration and creates the processing outputs.

The PP WS and the Reporter cooperate in order to save the produced outputs and to trace the Co-Registration procedures.

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SEGMENTATION

Descrizione


The PP WS gets the suitable inputs, performs the segmentation and creates the processing outputs.

The PP WS and reporter cooperate in order to save the produced outputs and to trace the segmentation procedure.

Diagramma

Created/modified by: Giorgio Cangioli
Creation Date: 10/09/2002

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PVE CORRECTION

Descrizione

The PP WS gets the suitable inputs, performs the PVE correction and generates the PVE corrected PET images and object values.


The PP WS and reporter cooperate in order to save the produced outputs and to trace the PVE correction procedure.

Diagramma

Created/modified by: **Giorgio Cangioli**
 Creation Date: **10/09/2002**
 Last Modified: **12/09/2002**
 Version: **0.2**

The PP WS allows users to:
 (1) select the algorithm to be applied.
 (2) choose the inputs images [TBD]
 (3) define the additional parameters needed for the selected algorithm [TBD].

To Be Defined / Verified :
 (1) Are the input images user selected ?
 (2) Which input parameters should be defined ?
 (3) Does this package use the segmentation and co-registration outputs as input images ?

The PP WS creates PVE corrected PET images/object values according to the selected algorithm.
 4 different techniques may be used for this purpose

The PP WS saves the PEV corrected outputs :
 (1) PVEC PET images
 (2) PVEC object values

TBD :
 (a) Should the outputs be saved on the Reporter together with the "tracing"report ?
 (b) How to save (file format, SOP Class) the PET images
 (c) How to save the numeric outputs (file format, SOP Class,...)

The PP WS makes use of the Reporter capabilities for tracing and recording the Partial Volume Correction procedure applied by using DICOM SR.

The SR shall include:
 (1) The algorithm selected
 (2) The input parameters used [TBD]
 (3) The input images (by reference)
 (4) [TBD how to (and if) manage analyze input/output]
 (5) The output images (by reference)
 (6) The output object values [TBD]

Even if the PVEC procedure may fail there is no interest in recording the reasons for failure.

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